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Beam	Treatment Un	Modality	Algorithm	Energy	Isocenter Location	X(cm)	Y(cm)	Z(cm)
1					Interest Point 1: 10cm	0.00	0.00	2.00

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The screenshot displays the Elekta treatment planning software interface. The top section shows a 2D cross-sectional view of a patient's head and neck with a treatment field defined by a red box and yellow lines. The bottom section is the 'Prescription' tab, which includes a table for defining the prescription. The table has three columns: 'Rx Dose (cGy)', 'Number of Fractions', and 'Fractional Dose (cGy)'. The values '100.0', '1', and '100.0' are entered in these columns, respectively. A red box highlights these three columns. Below the table, there are fields for 'Rescale' (set to 100.0 cGy) and 'Weight beams by' (set to Dose). At the bottom, a navigation bar shows 'Prescription' as the active tab, also highlighted with a red box.

Rx ID	Rx Site	Prescribe To	Rx Dose (cGy)	Number of Fractions	Fractional Dose (cGy)
Physician's Intent	A	Plan Isocenter	100.0	1	100.0

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The screenshot displays the Elekta treatment planning software interface. At the top right, a blue box highlights the text "10x10照射野の場合". The main workspace is divided into several panels: a top-left panel showing a patient's head and neck in a sagittal view with a red box indicating the target area; a top-right panel showing a grid of dose distribution; a bottom-left panel showing a cross-sectional view of the patient's head and neck with a red box indicating the target area; and a bottom-right panel showing a cross-sectional view of the patient's head and neck with a red box indicating the target area. The "Beams" tab is selected in the bottom left, and the "Geometry" sub-tab is selected in the bottom right. A table below the "Beams" tab shows the beam parameters for a 10x10 cm field.

Beam	Description	SSD (cm)	Gantry (deg)	Collimator (deg)	Couch (deg)	Asym	Width1 (cm)	Width2 (cm)	Length1 (cm)	Length2 (cm)
1	AP	90.01	0.0	0.0	0.0	2)	X1 5.00	X2 5.00	Y2 5.00	Y1 5.00

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