
ACCURAY

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
For Precision  
(Radixact System ver. 1.x)

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Create Machine QAを選択

Machine Name: 4010017 HRC\_6MV X-ray, 3.5MV MVCT  
Machine Serial Number: 4010017  
Machine Status: Ready  
Water Temperature: 38.90°C  
Service Mode: OFF  
Warmup Status: Complete

What's Next?  
Treatment Delivery Console  
Select one of the workflows below.

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
Thursday, 05 Apr 2018 11:56:56

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Basic Createを選択

The screenshot displays the Radixact ACCURAY software interface. At the top, the machine name is '4010017\_HRC\_6MV X-ray, 3.5MV, MVCT'. A 'What's Next?' box prompts the user to 'Create Machine QA' and select an option: 'create, view, or delete a Machine QA Procedure'. The 'Procedure' menu is open, showing several options: 'View', 'Basic Create', 'Create', 'Create From Patient Plan', and 'Delete'. The 'Basic Create' option is highlighted with a red box, and a red arrow points to it from a text box below. The text box contains the instruction: '“Basic Create”をクリックします。' (Click 'Basic Create').

Machine Name: 4010017\_HRC\_6MV X-ray, 3.5MV, MVCT

What's Next?  
Create Machine QA  
Select one of the options below to:  
create, view, or delete a Machine QA Procedure.

Print Tools Lock

User: training, training  
Version: 1.0.0.6

Procedure

View Query the Data Server or load a Procedure to view all its parameters.

**Basic Create** Define the basic parameters to create a Procedure.

Create View all parameters to define one or more as required to create a Procedure.

Create From Patient Plan View all parameters as taken from a patient plan to redefine one or more as required to create a Procedure.

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Radixact ACCURAY

Mode: Create Machine QA

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The screenshot shows a software interface for creating a procedure. At the top, there are buttons for "Save to XML File", "Save to Data Server", and "Cancel". Below these are several input fields and a list of procedures. The "Machine QA Name" field is set to "Output 1Gy". The "Procedures" list contains one entry: "5x10cm 10.3sec", with "Clone" and "Delete" buttons next to it. The "Description" field is also set to "5x10cm 10.3sec". The "Setup To Ready Longitudinal Distance (mm)" field is set to "0". The "Compression Factor" field is set to "10". The "Compression Type" dropdown menu is set to "BY\_COMPRESSION\_FACTOR". The "Beam" dropdown menu is set to "Treatment, Nominal MU: 1180, Pulse Rate: 300". At the bottom left, there is a "Fragment Summary" box with the following text: "Type of delivery: Static", "Total number of projections: 30", "Projections per second: 1.00", "Gantry Rotations: None", "Projections per rotation: 6.00", "Total Couch Distance: 0.00 mm", "Procedure Time: 30.00 seconds". Below the summary is a "Status" section which is currently empty.

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**Jaws**

Jaw Mode: Fixed

Front Jaw Position (mm): -21.0 Front Jaw: -21.0 Back Jaw: 21.0

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Procedure

Save to XML File Save to Data Server Cancel

Save the Procedure parameters to an XML file.

Machine QA Name: Output 1Gy

Procedures: 5x10cm 10.3sec

Description: 5x10cm 10.3sec

Setup To Ready Longitudinal Distance (mm): 0

Compression Factor: 10

Compression Type: BY\_COMPRESSION\_FACTOR

Beam: Treatment, Nominal MU: 1180, Pulse Rate: 300

Fragment Summary

Type of delivery: Static  
Total number of projections: 203  
Projections per second: 10.00  
Gantry Rotations: None  
Projections per rotation: 60.00  
Total Couch Distance: 0.00 mm  
Procedure Time: 20.30 seconds

Status

Look In: calibrationData

5x10 60sec CouchImages CTs JawProfiles sinograms StaticDDTSinograms TreatmentProcedures

Save Binaries

Checking this option will save copies of all binary data in the data definition to disk with the xml file

File Name: Output 1Gy

Files of Type: Accuray calibration files (\*.xml)

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
The screenshot shows the 'Procedure' configuration window. At the top, there are three buttons: 'Save to XML File', 'Save to Data Server', and 'Cancel'. A red box highlights the 'Save to Data Server' button, with an arrow pointing to it from a text box that says 'Save to Data Serverをクリックします。' (Click Save to Data Server).

Below the buttons, the 'Machine QA Name' is 'Output 1Gy'. The 'Procedures' list contains '5x10cm 10.3sec'. The 'Description' is also '5x10cm 10.3sec'. Other fields include 'Setup To Ready Longitudinal Distance (mm): 0', 'Compression Factor: 10', 'Compression Type: BY\_COMPRESSION\_FACTOR', and 'Beam: Treatment, Nominal MU: 1180, Pulse Rate: 300'.

A 'Fragment Summary' box shows: Type of delivery: Static, Total number of projections: 203, Projections per second: 10.00, Gantry Rotations: None, Projections per rotation: 60.00, Total Couch Distance: 0.00 mm, Procedure Time: 20.30 seconds.

A 'Status' section is visible at the bottom left.

Overlaid on the right is a dialog box titled 'Saving Machine QA Procedure to Data Server' with ID: 23208413. It lists the following steps: Validating parameters, Validating reference image file, Validating procedures, Validating Procedure: 5x10cm 10.3sec, Validating fragment [1] settings, Checking for matching beam, Validating fixed jaw settings, Validating sinogram file, Finished validating Procedure: 5x10cm 10.3sec, Finished validating parameters, Uploading to server..., Creating machine QA record, Building Machine QA, Loading Reference Image, Building procedure sets, Building procedure, Building fragment, Loading Sinogram, and Uploading files. At the bottom of the dialog, a message says 'Machine QA record created successfully' and a 'Close' button is highlighted with a red box. An arrow points from this box to a text box at the bottom of the slide that says 'Machine QA record created successfullyと表示されたらCloseをクリックします。' (Click Close when Machine QA record created successfully is displayed).

**Thank you**

ACCURAY