

ACCURAY

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For Planning Station  
(Hi-Art System, TomoTherapy H Series)

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Calibrate > Define > Procedure > Basic Create

The screenshot shows the TomoTherapy Operator Station interface. At the top, the patient information section displays "Patient: No Patient Selected" and various fields for DOB, Sex, ID, Oncologist, Disease, Plan Label, Plan Status, Plan Date, Position, and DQA Plan. Below this is a navigation bar with tabs for Scan, Register, Treat, Plan, and Calibrate. The Calibrate tab is selected. Underneath, there are sub-tabs for Define, Run, Procedure, Calibration, AOM, JAM, and Machine. The Define tab is selected, and within it, the Procedure sub-tab is active. The "Basic Create" button is highlighted with a red box. A red arrow points from the Calibrate tab to the Basic Create button. The interface also includes a "What's Next" help section, a "User Name: System User" field, and a "Current Machine: DiamondStandard" status.

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The screenshot shows the 'Define' tab of a procedure creation software. The interface includes several sections: 'Procedures' with a list containing '5x10cm', 'Number of Procedures' set to 10, 'Setup To Ready Longitudinal Distance (cm)' set to 0, 'Compression Factor' set to 10, and 'Compression Type' set to 'BY\_COMPRESSION\_FACTOR'. The 'AOM' dropdown is set to 'TX\_1CM\_BEAM, MU: 1180, TriggerRate: 300'. The 'Fragment Summary' section shows: Type of delivery: Static, Total number of projections: 221, Projections per second: 10.00, Gantry Rotations: None, Projections per rotation: 60.00, Total Couch Distance: 0.00 cm, Procedure Time: 22.10 seconds. The 'Status' section is empty. The right side of the interface shows 'Jaws', 'Gantry', 'Sinogram', and 'Couch' settings.

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Jaws、Gantry、Sinogramの設定

The screenshot shows a software interface for creating a procedure. The interface is divided into several sections: Scan, Register, Treat, Plan, and Calibrate. The 'Calibrate' section is active, and the 'Machine' tab is selected. The 'Current Machine' is 'DiamondStandard'. The 'Procedures' list is empty. The 'Jaws' section has 'Jaw Mode' set to 'Fixed' and 'Front Jaw: 2.1 Back Jaw: -2.1' selected. The 'Gantry' section has 'Gantry Mode' set to 'Fixed' and 'Gantry Start Angle (degrees)' set to '0.0'. The 'Sinogram' section has 'Sinogram Mode' set to 'Dynamic' and 'Projections Per Second' set to '10'. The 'Couch' section has 'Couch Mode' set to 'Fixed'. The 'Edit' button is highlighted.

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Current Machine: DiamondStandard

Define Run

Procedure Calibration AOM JAM Machine

Create the Procedure.

Save Save the Procedure parameters to an XML file.

Create Create the Procedure on the Data Server.

Cancel Cancel to choose a different task.

Patient Name: Output

Procedures: 5x10cm Clone Delete

Description: 5x10cm

Number of Procedures: 10

Setup To Ready Longitudinal Distance (cm): 0

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The screenshot shows a software interface for creating a procedure. The main window has tabs for 'Scan', 'Register', 'Treat', 'Plan', and 'Calibrate'. Under 'Calibrate', there are sub-tabs for 'Define' and 'Run'. The 'Define' sub-tab is active, showing a 'Procedure' section with a 'Save' button highlighted by a red box and an arrow pointing to it. A red callout box contains the text: "Save"をクリックします。

Below the 'Save' button, there are fields for 'Patient Name' (Output), 'Procedures' (5x10cm), 'Description' (5x10cm), 'Number of Procedures' (10), 'Setup To Ready Longitudinal Distance (cm)' (0), 'Compression Factor' (10), 'Compression Type' (BY\_COMPRESSION\_FACTOR), and 'AOM' (TX\_1CM\_BEAM, MU: 1180, TriggerRate: 300).

A 'Save' dialog box is open in the foreground. It shows the 'Save In' location as 'calibrationData' and a list of folders: 'CTs', 'sinograms', 'StaticDDTSinograms', and 'TreatmentProcedures'. The 'File Name' field contains '5x10cm 1Gy'. The 'Files of Type' dropdown is set to 'TomoTherapy calibration files (\*.xml)'. A red callout box points to the 'File Name' field with the text: File Nameを入力します。 ※ 本例では“Output 1 Gy”です。

The 'Save' dialog box has a 'Save' button highlighted by a red box and an arrow pointing to it. A red callout box contains the text: "Save"をクリックします。

In the bottom left corner, there is a 'Fragment Summary' section with the following details:

- Type of delivery: Static
- Total number of projections: 221
- Projections per second: 10.00
- Gantry Rotations: None
- Projections per rotation: 60.00
- Total Couch Distance: 0.00 cm
- Procedure Time: 22.10 seconds

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Current Machine: DiamondStandard

Save Save the Procedure parameters to an XML file.

**Create** Create the Procedure on the fly.

Cancel Cancel to choose a different task.

Patient Name: Output

Procedures: 5x10cm

Description: 5x10cm

Number of Procedures: 10

Setup To Ready Longitudinal Distance (cm): 0

Compression Factor: 10

Compression Type: BY\_COMPRESSION\_FACTOR

AOM: TX\_1CM\_BEAM, MU: 1180, TriggerRate

**Fragment Summary**

Type of delivery: Static  
Total number of projections: 221  
Projections per second: 10.00  
Gantry Rotations: None  
Projections per rotation: 60.00  
Total Couch Distance: 0.00 cm  
Procedure Time: 22.10 seconds

**Status**

Validating parameters  
Converting data definition model to XML  
Validating XML schema  
Validating Data Definitions  
Validating binary file: PlanningCT/TomoImage/Filename  
Validating binary file: PlanningCT/KVCT\_Image/Filename  
Validating procedure sets  
Validating Procedure: 5x10cm 60sec  
Checking for matching AOM  
Validating scheduled procedure [1] settings  
Validating binary file: mlcProcedure/Filename  
Validating fixed jaw settings  
Checking for matching JAM  
Finished validating Procedure: 5x10cm 60sec  
Checking for empty filename nodes  
Checking for current calibration  
Finished validating parameters

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